

Edwards Pirani 501 Vacuum Gauge with PRE10K Sensor



Range: Atmospheric to 0.001 mbar

Price: ~~Not Listed~~ \$177.50

Discount \$-177.50

[Ask a question about this product](#)

Manufacturer [Edwards](#)

Edwards Pirani 501 Vacuum Gauge with PRE10K Sensor:

- Features analog display
- Connection flange: KF10
- Range: Atmosphere to 0.001 mbar
- Includes black, 6 ft. AC power cord

**ANALYTICAL
INSTRUMENTS**

Sales - Service - Buy

Warrantied Lab Equipment

**ANALYTICAL
INSTRUMENTS
LLC**

**Sales | Service | Buy
Refurbished Laboratory Equipment
With Warranty**

- Operates on 120 VAC, 50\60 Hz, single phase main power
(Or 200\230 VAC, 50\60 Hz, single phase main power)
-
