


Eppendorf 5417C Tabletop Centrifuge



Maximum Speed: 14,000 RPM

Rating: Not Rated Yet

Price

Sales price \$2,200.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Eppendorf 5417C Tabletop Centrifuge:

- Automatic rotor detection
- Max. Capacity: 30 x 1.5 – 2.0ml
- Maximum Speed: 16,400 RPM (fixed rotor)
- Maximum Speed: 10,000 RPM (swinging bucket)
- Digital LED readout and display
- Display allows users to select time & speed, in variable seconds or minutes & RPM or RCF units
- Includes black 6' main AC power cord
- Simple easy use keypad & push-button operation
- Compatible rotors available for purchase (not included); Options can include Eppendorf 30-place or 24-place rotors, swinging bucket 8-place, or six 8-tube strips
- Small notch of missing plastic on the front left of chassis
- Operates on 120 VAC, 60 Hz, single phase main power
- Includes 90-day warranty directly from A.I.