

International Cryogenics Model IC5D Dewar Liquid Nitrogen Storage



5.0 Liter storage LN2 dewar

Rating: Not Rated Yet

Price

Sales price \$295.00

Tax amount

[Ask a question about this product](#)

Description

International Cryogenics Model IC5D Dewar Liquid Nitrogen Storage:

- Double-wall construction with thermal barrier, insulation and vacuum
- Storage dewar for Liquid Nitrogen (LN2)
- Outside Diameter: 8.75" (222 mm)
- Neck Diameter (ID): 2.25" (57 mm)
- Total Height: 15.75" (400 mm)
- Capacity: 5.0 Liters
- Full weight: 15.0 lbs (7.0 kg.)
- Empty weight: 6.0 lbs (3.0 kg.)
- Static hold time: 25 days
- Evaporation rate: 0.2 Liters \ day