

MVE Tech 50 Liquid Nitrogen LN2 Tank



Capacity 50 Liters

Rating: Not Rated Yet

Price

Sales price \$2,485.00

Tax amount

[Ask a question about this product](#)

Manufacturer [MVE](#)

Description

MVE Tech 50 Liquid Nitrogen LN2 Tank:

- Capacity 50 Liters
- Uses two pressure gauges
- Includes vent valve to prevent pressure build-up
- Includes a fill port and an electric controlled outlet port
- Electric controlled outlet port uses 120 VAC main power
- Includes 90-day warranty directly from A.I.