

Taylor Wharton 35VHC Liquid Nitrogen Dewar



Capacity: 35 Liter

Rating: Not Rated Yet

Price

Sales price \$380.00

Tax amount

[Ask a question about this product](#)

Manufacturer [Taylor Wharton](#)

Description

Taylor Wharton 35VHC Liquid Nitrogen Dewar:

- Capacity: 35 Liters
- Includes roller base
- Holds six (6) canisters
- Strudy design, high capacity
- LN2 static hold time of 130 days
- Evaporation rate: 0.27 Liters per day